

# Bacti-cinerator Sterilizer

## IS800-A



### Description:

- ◆ IS800-A Bacti-cinerator Sterilizer is used in degassing for inoculating loops and needles, instead of alcohol lamps. Due to its convenience, sterilizer is widely used in the environment of biosafety cabinets, clean benches. They are also used besides exhaust fans or on mobile vans.

#### Application:

Widely used in biosafety cabinets, clean bench, near exhaust fan, mobile vehicles, even outdoor and other other severe environment can be easily used.

### Features:

- ◆ Safety, convenient. Totally take the place of alcohol lamps.
- ◆ Only 5 seconds to 7 seconds are needed, because of its highest temperature of 825°C.
- ◆ Smart, light, easy to handle.
- ◆ Used in anaerobic chamber is allowed.
- ◆ Podzolic organics prevent pollutions effectively.
- ◆ Smart technique makes it possible to use for a long time and much safer.
- ◆ 12-month warranty.

#### Parameters:

Project	Parameters
Max. Temp. in central zone	825°C±25°C
Length of heating zone	140mm
Ambient temperature	5°C~40°C
Voltage	AC220V~ 50-60Hz
Dimension(mm)	120×155×180

Project	Parameters
Max. Out diameter of drug elimination	Φ15mm
Heating time	15min
Fuse	250V 3A 5×20
Power	150W
Weight(kg)	1.0



Hangzhou Rui Cheng Instrument CO., LTD

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China  
TEL: 86-571-88630940 FAX: 86-571-85126370 Email:sales2@ruiwell.com  
<http://www.ruiwell.net>