

Biological Indicator Incubator

B80



Description:

- ◆ The B80 Biological Indicator Incubator is compact, versatile, reliable and easy to use. It is designed for single temperature incubation of the EO Biological Indicators. The thermometer is included with the incubator, providing easy verification of incubation temperature.

Features:

- ◆ Energy-efficient design.
- ◆ Applicable to 24V power adapter.
- ◆ Accurate to $\pm 0.3^{\circ}\text{C}$.
- ◆ Process of experiment is visible.
- ◆ Electronic calibration. User can reset to any local standard.
- ◆ Power supply provided is UL, CSA and CE approved.
- ◆ Custom blocks made to fit any sample holder.
- ◆ 12-month warranty.

Parameters:

Project	Parameters
Temp. control range	RT+5°C~100°C
Timing range	1min~99h
Display accuracy	0.1°C
Heating time	≤6min(20°C~56°C)
Voltage	DC24V/1.67A

Project	Parameters
Temp. accuracy	≤±0.3°C
Temp. uniformity	≤±0.3°C
Temp. stability	±0.3°C
Dimension(mm)	145×113×110
Weight(kg)	0.5

Optional Blocks:

Type	Diameter	Capacity
MD40	Φ 10mm	15
MD41	Φ 12mm	12
MD42	Φ 15mm	8

Type	Diameter	Capacity
MD43	Φ 16mm	8
MD45	Φ 9mm	18
---	---	---



Hangzhou Ruicheng Instrument Co., Ltd

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China
TEL: 86-571-88630940 FAX: 86-571-85126370 Email: sales2@ruiwell.com
<http://www.ruiwell.net>