


# Shaker For Microplates

## MS100



### Description:

- ◆ MS100 Shaker for Microplates is a type of high speed shaker, which is with brushless DC motor and microprocessor controlled technique, specially for elisa plates(96/384 wells), tissue culture plates(24/ 48/96 wells).

### Features:

- ◆ LCD display system status and parameters.
- ◆ Stable and reliable shaking.
- ◆ Easy to operate with one touch knob.
- ◆ Automatic fault detection and buzzer alarm function.
- ◆ 12-month warranty.

### Parameters:

Project	Parameters
Time range	1min ~ 99h59min
Mixing orbit	3mm
Sample capacity	4plates
Voltage	AC110V~240V 50/60 Hz
Weight(kg)	4.0

Project	Parameters
Speed range	200~1600 rpm
Speed step value	10rpm
Power	45W
Dimension(mm)	280×270×140
---	---



Hangzhou RuiCheng Instrument CO., LTD

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China  
TEL: 86-571-88630940 FAX: 86-571-85126370 Email:sales2@ruiwell.com  
<http://www.ruiwell.net>