

Mini Vortex Mixer

MIX1000



Description:

- ◆ MIX1000 Mini Vortex Mixer is designed for “touch” mixing a variety of test tubes. With a 4.5mm orbit and adjustable speed of 1000-2800rpm, the mini vortex mixer instantly vortexes even the largest samples, including nearly full 50ml tubes. The small footprint allows it to be used on even the most crowded bench.

Features:

- ◆ Shaking speed: 1000-3500rpm.
- ◆ Powerful and reliable motor designed for daily laboratory use.
- ◆ “Touch” pressure activated operation.
- ◆ Accept a variety of test tubes, including the tubes up to 50ml.
- ◆ Mini-sized, allowing it to fit on even the most crowded bench.

Parameters:

Project	Parameters
Speed	1000-3500rpm
Type of movement	Orbital
Max test tube size	Φ30mm
Voltage	AC100-240V, 50/60Hz
Dimension(mm)	125×115×76

Project	Parameters
Type of motor activation	Touch pressure activated operation
Orbit	4.5mm
Max load	0.1kg
Power	24W
Weight(kg)	0.55



Hangzhou Ruicheng Instrument Co., Ltd

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China
TEL: 86-571-88630940 FAX: 86-571-85126370 Email: sales2@ruiwell.com
<http://www.ruiwell.net>