

# Biological Indicator Incubator

## B80



### Description:

- ◆ The B80 Biological Indicator Incubator is compact, versatile, reliable and easy to use. It is designed for single temperature incubation of Steam or EO Biological Indicators.

### Features:

- ◆ Energy-efficient design.
- ◆ Applicable to 24V power adapter.
- ◆ Accurate to  $\pm 0.3^{\circ}\text{C}$ .
- ◆ Process of experiment is visible.
- ◆ Electronic calibration. User can reset to any local standard.
- ◆ Power supply provided is UL, CSA and CE approved.
- ◆ Custom blocks made to fit any sample holder.

### Parameters:

Project	Parameters
Temp. control range	RT+5°C~100°C
Timing range	1min~99h
Display accuracy	0.1°C
Heating time	≤6min( 20°C~56°C )
Voltage	DC24V/1.67A

Project	Parameters
Temp. accuracy	≤±0.3°C
Temp. uniformity	≤±0.3°C
Temp. stability	±0.3°C
Dimension(mm)	145×113×110
Weight(kg)	0.5

### Optional Blocks:

Type	Aperture	Number of test tube
MD40	Φ 10mm	15
MD41	Φ 12mm	12
MD42	Φ 15mm	6

Type	Aperture	Number of test tube
MD43	Φ 16mm	6
MD45	Φ 12mm(cone-shape bottom)	14
MD45	Φ 9mm	18



Hangzhou Rui Cheng Instrument CO., LTD

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China  
TEL: 86-571-88630940 FAX: 86-571-85126370 Email:sales2@ruiwell.com  
<http://www.ruiwell.net>