

# Mini Vortex Mixer

## MIX1000



### Description:

- ◆ MIX1000 Mini Vortex Mixer is designed for “touch” mixing a variety of test tubes. With a 4.5mm orbit and adjustable speed of 1000-3500rpm, the mini vortex mixer instantly vortexes even the largest samples, including nearly full 50ml tubes. The small footprint allows it to be used on even the most crowded bench.

### Features:

- ◆ Shaking speed: 1000-3500rpm.
- ◆ Powerful and reliable motor designed for daily laboratory use.
- ◆ “Touch” pressure activated operation.
- ◆ Accept a variety of test tubes, including the tubes up to 50ml.
- ◆ Mini-sized, allowing it to fit on even the most crowded bench.

### Parameters:

| Project            | Parameters          |
|--------------------|---------------------|
| Speed              | 1000-3500rpm        |
| Type of movement   | Orbital             |
| Max test tube size | Φ30mm               |
| Voltage            | AC100-240V, 50/60Hz |
| Dimension(mm)      | 125×115×76          |

| Project                  | Parameters                         |
|--------------------------|------------------------------------|
| Type of motor activation | Touch pressure activated operation |
| Orbit                    | 4.5mm                              |
| Max load                 | 0.1kg                              |
| Power                    | 24W                                |
| Weight(kg)               | 0.55                               |



Hangzhou Ruicheng Instrument Co., Ltd

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China  
TEL: 86-571-88630940 FAX: 86-571-85126370 Email: sales2@ruiwell.com  
<http://www.ruiwell.net>