

Orbital Shaker

MPS-10



Description:

- ◆ MPS-10 Orbital Shaker is a high speed shaker with direct brushless DC motor and microprocessor controlled technique. It can be applicable in different laboratories: in microbiology, biotechnology, medical analysis, etc.

Features:

- ◆ LCD displays time and shaking speed.
- ◆ Easy to operate with one-touch knob.
- ◆ Warm root and variety of platforms can be chosen.
- ◆ When time control ends, instrument will make alarm voice.
- ◆ Gentle and reliable shaking with direct brushless DC motor.

Parameters:

Project	Parameters
Time Range	1min ~ 99h59min
Amplitude	10mm
Power	50W
Dimension(mm)	280×270×110

Project	Parameters
Speed Range	50~300rpm
Max. load	2.5kg
Voltage	AC110V~240V 50/60Hz
Weight(kg)	7.0

Optional Platforms:

Type	Description
PP-1	Flat platform for 8pcs petri dish
PP-2	Universal platform with adjustable bars, for 6x100ml or 4x250ml conical flask

Type	Description
PP-3	Platform with clamps for 12x100ml conical flask
PP-3	Platform with clamps for 6x250ml conical flask



Hangzhou Rui Cheng Instrument CO., LTD

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China
TEL: 86-571-88630940 FAX: 86-571-85126370 Email:sales2@ruiwell.com
<http://www.ruiwell.net>